Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

GC Fuji TEMP LT

GLASS IONOMER PROVISIONAL LUTING CEMENT

Prior to use, carefully read the instructions for use.

EN

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Temporary (provisional) cementation of all types of all-ceramic, resin and metal-based crowns and bridges
2. Temporary (provisional) cementation of all types of all-ceramic, resin and metal-based crowns and bridges on implant abutments

CONTRAINDICATIONS

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

IDENTIFICATION OF PARTS (Fig. 1)

A. Paste Pak Cartridge

1. Cover 2. Inner plunger 3. Cartridge tip

B. Paste Pak Dispenser

4. Piston 5. Adjusting guide 6. Lever 7. Sliding block 8. Piston release

PASTE PAK CARTRIDGE LOADING

1. Make sure that the piston is completely retracted into the dispenser (Fig. 2).

NOTE:

To retract the piston, push the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.

2. Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the mark (▲) on the cartridge lines up with the mark (▲) on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.

3. After cartridge loading, push the piston release forward until it stops (Fig. 4).

- a) Push the piston release forward.
- b) Do not push the upper side.

NOTE:

To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.

DIRECTIONS FOR USE

Paste / Paste ratio (g/g)	0.7 / 1.0
Mixing time (sec.)	10"-15"
Working time (min., sec.) (at 23°C, from start of mixing)	2'00"
Final finishing commencing time (min., sec.) (after removing excess)	3'00"

1. ABUTMENT PREPARATION

- 1) Clean the tooth/abutment surfaces with pumice and water. For pulp capping, use calcium hydroxide.
- 2) Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should be moist (glistening). Overdrying may cause difficulties in cleaning the residual cement after debonding.

NOTE:

- a. When a moderate retention is desired and especially when the fit of restoration is tight, a thin layer of separation medium like petroleum jelly can be applied on the abutment surface. Prior to final, permanent cementation remove petroleum jelly with alcohol first then water, as petroleum jelly may hinder setting or bonding of the permanent cement.
- b. The use of desensitizer may lower the bonding of GC Fuji TEMP LT.

2. RESTORATION PREPARATION

It is recommended to sandblast the internal surface of the restoration when increased retention is desired. For sandblasting procedure, refer to the respective manufacturer's instructions. Before cementation, apply a thin layer of petroleum jelly to the outer marginal surface of the restoration. This will ease the clean-up of excess cement.

3. DISPENSING

- 1) Removing the cartridge cover.
- 2) Depress the lever to dispense the required amounts of paste onto the mixing pad.

NOTE:

- a. Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge. The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 5).

1. Position of adjusting guide for one inlay.
2. Position of adjusting guide for one crown.

- b. The paste may be dispensed in the shape of a rod, but this does not affect the quality of the paste.

- 3) Depress the lever while moving the cartridge and dispenser backward (Fig. 6). Keep depressing the lever (Fig. 7-1), level and cut off the extruded material from the cartridge tips by moving the cartridge and dispenser to an upright position on the mixing pad (Fig. 7-2).

- 4) When the lever is released after dispensing, residual pastes at the cartridge tips are retracted into the cartridge. Clean the cartridge tip and replace the cartridge cover.

NOTE:

The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser.

4. MIXING

Spread the two pastes thinly on the mixing pad using a plastic spatula. Mix vigorously and thoroughly with lapping strokes, for 10 seconds. Be careful not to incorporate air bubbles.

NOTE:

- a. Larger amount of paste requires a longer mixing time (15 seconds).
- b. Do not mix GC Fuji TEMP LT with other materials.

5. CEMENTATION

- 1) Coat the fitting surface with a thin layer of mixed material. The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time.
- 2) Seat the restoration.
- 3) To remove excess cement, allow the cement to self-cure until excess cement feels rubbery (approx. 1 minute).
- 4) Remove excess cement with a scaler, probe and/or floss.
- 5) Maintain isolation until the cement has set. Finishing can be started approx. 3 minutes after excess cement removal.

CLEANING / STERILIZATION OF PASTE PAK DISPENSER

For cleaning and sterilization of PASTE PAK DISPENSER, refer to the respective instructions for use.

SHADE

Universal

STORAGE

Store at room temperature away from heat, moisture and direct sunlight (4-25°C)(39.2~77.0°F) (Shelf life: 2 years from date of manufacture)

PACKAGE

- 1) GC Fuji TEMP LT Cartridge (13.3g / 7.2mL) (2), mixing pad (1)
- 2) Paste Pak Dispenser (1)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and / or cocoa butter can be used to isolate the operation field from the skin or oral tissue. Wear plastic or rubber gloves for protection.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Avoid getting material on clothing.

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