Prior to use, carefully read
the instructions for use.

EN

PERI COMPOUND

BORDER MOULDING IMPRESSION MATERIAL

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

Making preliminary impressions (by the border moulding technique for dentures).

DIRECTIONS FOR USE

1. Soften material using a Bunsen burner or an alcohol lamp. Apply to the periphery of the custom tray (Fig. 1).
2. Immerse the tray in 45°C water to obtain plasticity of the material (Fig. 2).
3. Insert in the mouth, let the patient have sucking and swallowing motions with the help of the finger in order to perform the border moulding (Fig. 3).
4. Remove any excess material using a scalpel to complete the border moulding (Fig. 4).
5. For the next steps of the procedure, follow standard techniques.

STORAGE

Recommended for optimal performance, store at normal room temperature (4-25°C / 39.2-77.0°F).

PACKAGE

Box of 125 g (15 sticks)

NOTE

1. When softening material, be careful not to get burned. Make sure that the softened material has an appropriate temperature before insertion of the tray in the patient's mouth.
2. The obtained impression should be cleaned, then disinfected using an appropriate disinfectant according to manufacturer's instructions.
3. Do not soften material at an excessively high temperature to avoid stickiness of material.
4. When it is difficult to remove excess of softened material from the denture surface, leave alone for a while until it hardens. Then remove with a scalpel and clean with benzene.
5. Do not mix with other compound materials.
6. In rare cases this product may cause sensitivity to some persons. If such reactions are experienced, discontinue use of the product and refer to a physician.
7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<http://sea.gcasiadental.com>
<http://www.gcaustralasia.com>

They can also be obtained from your supplier.

Last revised: 12/2019

GC

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza, Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

MANUFACTURED by
NIMI CHEMICAL INDUSTRIAL CO., LTD.
4-13 Yokoyama-cho, Kiryu, Gunma 378-0055, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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